

Service Dog Designs

by Bold Lead Designs, llc

Walking Assistance/ Basic Mobility Harness

\$285 for average dogs (60-100 pounds)
\$315 for extra large dogs (over 100 pounds)

A complete harness includes:

- Harness saddle with straps custom made to fit your dog
- Fixed leather handle made to your desired height
- Removable and adjustable pull-strap
- 7 D-rings on the back and sides of harness
- Great Fit Guarantee: fit adjustment if necessary within 4 weeks of delivery
- Lifetime customer support

Modifications are available; just tell us what you need.

Extras:

ID Patches & Badges: various designs and phrases available

Double Patches –2-sided badges Velcro onto leash or harness strap \$16

Clip on buttons—various designs available \$6

Harness Pouch: a nifty accessory pouch to hold essentials, \$24
like ID, money, phone, treats, medicine, cleanup bags, etc.

We offer a variety of leashes, collars, and training accessories.

Check out Bold Lead Designs' main website to see all that we offer:
www.boldleaddesigns.com

Replacement Parts:

Our harnesses have been designed with interchangeable parts for the greatest flexibility. Perfect for adapting a harness to fit a new dog.

Breast plate with sheepskin lining (includes all straps to attach to harness): \$40

Belly strap (add \$10 for sheepskin lining): \$20

Girth (with sheepskin lined pad): \$20

Pull-Strap: \$25

Shipping:

We ship anywhere! Standard shipping is insured USPS Priority Mail (\$15-18 for a complete harness). Ask for a quote for other shipping options and international orders (outside the US).

If you are interested in a harness with a rigid brace handle, please see the information for our Mobility Support Harness.

Order Form:

You may order by phone, snail mail, e-mail, or online.

We accept payment by check, PayPal, or credit card.

Make payment to Bold Lead Designs.

When ordering, please provide the following information:

Date:

Contact Name:

Institution/Trainer Name:

Mailing Address:

Phone Number:

Email:

Dog's Name:

Dog's Breed:

User's Name:

Items you would like to order:

Dog's Measurements:

- 1) Height _____
- 2) Girth _____
- 3) Belly _____
- 4) Width _____
- 5) Dog's weight _____ pounds

Fixed Handle height _____
(from dog's back to handle)

Notes on modifications or features:


www.servicedogdesigns.com

- **Katrina Boldry, Owner & Designer • email kboldry@boldleaddesigns.com •**
- **Phone: 303-856-3012 • Mailing address 393 Lima Street, Aurora, CO 80010 •**

Service Dog Designs is a division of Bold Lead Designs, LLC

How to measure your dog:

- Have the dog stand square and balanced. You may want a friend to help you.
- Measure **height** from withers to the floor (lay a book across the shoulders if you like).
- For girth and belly measurements, use a cloth measuring tape or a string. Make a full circle around the dog's body. Snug the tape or string up tight for accurate measurements (we want the dog's body size, underneath the fur).
- Locate the "sweet spot": this is the dip in the dog's back just behind the withers (the bony shoulder blades).
- Measure the **girth** from the "sweet spot" around the deepest part of the dog's chest.
- Measure the **belly** from the "sweet spot," around the soft part of the belly, at a diagonal. Be sure the tape is pulled up tightly **BEHIND** the dog's rib cage (flanks) and then back to the "sweet spot."
- Measure the **width** straight across the **FRONT** of the dog's chest. This is easiest with a regular flat ruler. Sit in front of your dog, and hold the ruler in both hands across the dog's breastbone. Slide your fingers in until you are touching each side of the dog's shoulders. (Imagine that your dog is squeezing through a tiny doorway; what is the narrowest space he can squeeze his shoulders through?)
- Make these measurements 2 or 3 times for accuracy. Each time may give you different results, don't panic; just average the numbers.
- Remember that the harness will be adjustable!



How to determine the balance handle height:

Please indicate the desired height for the balance handle. The balance handle is the leather loop that sticks straight up off the harness saddle. Remember that the fixed leather balance handle cannot be adjusted later, as it is permanently attached to the harness saddle.

Note: The removable pull-strap will be made adjustable and can be clipped in multiple positions on the harness.

Option #1:

- Stand and loosely hold a pencil in your hand
- Allow your arm to hang naturally with the elbow slightly bent, keep the arm in a position that will be comfortable for balancing against your dog.
- Measure from the pencil to the floor
- Subtract the dog's height to equal the handle measurement

Option #2: the best option

- Indicate a specific handle height requirement (most are around 4 or 6 inches, depending on your dog's height)
- Clearly indicate this measurement when ordering.

Confused? Have questions? Don't worry. Katrina is here to help.

Call 303-856-3012

Email kboldry@boldleaddesigns.com